# ***Integra RNA/DNA MAG extraction kit***

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## EXMB003 INTEGRA RNA/DNA EXTRACTION KIT

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

### 1. GENERAL INFORMATION

#### 1.1 Intended use

**GENESPARK INTEGRA RNA/DNA MAG EXTRACTION KIT** has been designed for rapid manual or automated extraction of nucleic acids from a wide range of samples (eg human samples, veterinary samples, insects, food samples, etc.). The kit is designed to be used with KingFisher™ Flex magnetic particle processor or other magnetic separation systems. Purified RNA / DNA can be used directly as a template for RT-PCR, PCR or any type of enzymatic reactions. This kit allows easy automation on common liquid treatment instruments or automatic magnetic separators. The actual procedure time depends on the configuration of the instrument and the magnetic separation system used.

#### 1.2 Principle

**GENESPARK INTEGRA RNA/DNA MAG EXTRACTION KIT** is designed for nucleic acid analysis by polymerase chain reaction (PCR) or RT-PCR, where analyte isolation from various sample materials is required. To this end, the sample is treated by incubation in a proprietary lysis buffer.

To bind nucleic acids to magnetic particles, the same INTEGRA particles are added to the lysate. In the presence of a salt, the nucleic acids bind specifically to the surface of the particles INTEGRA.

After magnetic separation, paramagnetic beads are washed to remove contaminants and salts using **Inhibitor Removal buffer** (P2) and the **Wash buffer** (P3). The residual ethanol from the previous washing steps is removed by airdrying. Finally, pure RNA / DNA is eluted with **Elution buffer** (P4) low in salt. Purified RNA / DNA can be used directly for applications as described below.

### 2. KIT CONTENTS

**INTEGRA RNA/DNA MAG EXTRACTION KIT** (REF EXMB003-50) can execute 50 extractions.

**INTEGRA RNA/DNA MAG EXTRACTION KIT** (REF EXMB003-200) can execute 200 extractions.

Tavola 1 – Kit contents

REF	EXMB003-50	EXMB003-200
<b>PV1 Binding buffer</b>	1 x 30 ml + 26 ml 2-propanol	2 x 30 ml + 26 ml 2-propanol
<b>PA Poly A</b>	2 mg	4 mg
<b>P2 Inhibitor Removal buffer</b>	1 x 33 ml + 20 ml absolute ethanol	2 x 33 ml + 20 ml absolute ethanol
<b>P3 Wash buffer</b>	2 x 10 ml + 40 ml absolute ethanol	2 x 20 ml + 80 ml absolute ethanol
<b>P4 Elution buffer</b>	1 x 10,5 ml	1 x 21 ml
<b>MB Magnetic beads</b>	2 x 1.0 ml	4 x 1.0 ml

All solutions are clear and should not be used when precipitation is formed. Heat solutions at +18 to +25 ° C or in a 37 ° C water bath until the precipitates have melted.

## EXMB003 INTEGRA RNA/DNA EXTRACTION KIT

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

### 3. STORAGE INSTRUCTIONS

All the components of **GENESPARK INTEGRA RNA/DNA MAG EXTRACTION KIT** must be stored at +18 to + 25 ° C. If they are stored correctly, all kit components are stable until the expiration date printed on the label. Please note that improper storage at +2 to + 8 ° C (refrigerator) or -15 to -25 ° C (freezer) will adversely affect the purification of the nucleic acid when precipitation is formed in the solutions. Therefore the INTEGRA kits are always shipped at +18 to + 25 ° C. The reconstituted Poly A RNA Work Solution should be aliquoted. The aliquots are stable for 12 months if stored at  $\leq -18$  ° C.

### 4. UNPROVIDED MATERIALS

- Use of **INTEGRA** con *KingFisher™ Flex Magnetic Particle Processor* e con *GENESPARK MD-EX Magnetic particles Processor (REF EQMD005)* it is recommended but the kit can also be used with other processors.
- Microtubes for centrifugation Nuclease-free 1.5 o 2.0 ml.
- Magnetic separation plate, for example Square-well Block (96-well block with 2.1 ml square-wells)
- Elution plate for collecting purified nucleic acids, for example Elution Plate V-bottom (96-well microtiterplate with 0.3 ml u-bottom wells)
- Pipettes with tips for molecular biology or Tip Comps (e.g. KingFisher 96 tip comb for magnetic separation).
- microcentrifuge bench-top 13,000 x RPM .
- Absolute ethanol
- Thermoblock or laboratory oven (only for DNA extraction)

Note: consumables for specific automatic processors are not included in the kit are dependent of the mode of sample preparation.

**For manual applications we suggest the adoption of the TABLET PRISMA TUTOR for the standardization of molecular biology methods in manual version.**

### 5. WARNING AND PRECAUTIONS

1. This Kit should only be used by qualified personnel.
2. Good Laboratory Practice (GLP) must be applied.
3. Clinical samples should always be considered as potentially infectious material and all equipment used must be treated as potentially contaminated.
4. The integrated release buffer and inhibitor removal buffer contain hydrochloride guanidine which is an irritant. Always wear gloves and follow standard safety precautions to minimize contact during handling. Do not use bleaches for waste treatment!
5. Do not let these swabs touch the skin, eyes or mucous membranes. In case of contact, immediately flush the affected area with large quantities of water; otherwise the reagent can cause burns. If the reagent is poured, dilute the spill with water before cleaning it.
6. Do not store or use tampons near human or animal food.
7. Always wear gloves and follow standard safety precautions when handling these swabs
8. Exercise the normal precautions necessary for the handling of all laboratory reagents.

**EXMB003 INTEGRA RNA/DNA EXTRACTION KIT**

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

9. Do not collect reagents from different lots or from different bottles of the same lot. Immediately after use, close all bottles to avoid leakage, changes in buffer concentrations or buffer conditions. After the first opening, keep all the bottles in an upright position.
10. Do not use a kit after the expiration date.
11. Avoid contact of the release buffer and inhibitor removal with skin, eyes or mucous membranes. In case of contact, wash immediately with large quantities of water. Possible burns can occur if not treated. If the reagent is leaked, dilute with water before drying it dry.
12. Do not use any modified ethanol.
13. Use only calibrated pipettes
14. All originating materials and all resulting waste should be considered potentially infectious. Clean and disinfect all work surfaces perfectly with disinfectants recommended by local authorities.
15. Do not eat, drink or smoke in the laboratory work area.
16. Do not pipette by mouth.
17. Wear protective disposable gloves, laboratory coats and eye protection when handling specimens and kit reagents.
18. Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent bottles.
19. Sterile sterile pipettes are recommended.
20. Wash hands thoroughly after taking samples and test reagents.
21. Dispose of unused reagents and waste should be performed in accordance with national, federal and local regulations.
22. Material Safety Data Sheets (MSDS) are available on request from Immunospark.

**6. OPERATING PROCEDURE**

**6.1 Reagents preparation**

The following table summarizes the conditions for the preparation of the various components of the procedure and the consequent storage conditions:

**Table 2 – Reagents preparation**

REF	Reconstitution/Preparation		Storage and stability	Scope
	EXMB003-50	EXMB003-200		
<b>Poly A (PA)</b>	Dissolve in 0.5 ml of Elution Buffer and prepare aliquots of 50 ul.	Dissolve in 1.00 ml of Elution Buffer and prepare aliquots da 50.	Store aliquots at $\leq -18^{\circ}$ C, stable for 12 months.	Additive of the binding buffer.
<b>Binding Buffer (PV1)</b>	Add 26 ml 2-propanol Binding Buffer and mix well. Label and date the bottle.	Add 26 ml 2-propanol Binding Buffer and mix well. Label and date the bottle.	Store at $+ 18^{\circ}$ C - $25^{\circ}$ C. stable until expiration.	Lysis of cells
<b>Inhibitor Removal Buffer (P2)</b>	Add 20 ml of absolute ethanol to the Residue Removal Buffer and mix well. Label and date the bottle	Add 20 ml of absolute ethanol to the Residue Removal Buffer and mix well. Label and date the bottle	Store at $+ 18^{\circ}$ C - $25^{\circ}$ C. stable until expiration.	Removal of PCR inhibitors from nucleic acids.

**EXMB003 INTEGRA RNA/DNA EXTRACTION KIT**

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

<b>Wash Buffer (P3)</b>	Add 40 ml of absolute ethanol to the Wash Buffer and mix well. Label and date the bottle	Add 80 ml of absolute ethanol to the Wash Buffer and mix well. Label and date the bottle	Store at + 18 ° C - 25 ° C. stable until expiration.	Removal of salts, proteins and other impurities from nucleic acids.
-------------------------	------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------	------------------------------------------------------	---------------------------------------------------------------------

DNA samples can be purified from a wide range of sample sources:

**Table 3 – possible sample sources**

Sample type	Volume/Quantity	Volume Binding Buffer	Pre-treatment of the sample
<b>Stool</b>	Pea-size	500 ul	Prepare a suspension in 1.5 ml of sterile dH2O. Vortex and filter the sediments. Use 200 µl of supernatant.
<b>Swab</b>		500 ul	
<b>Liquid samples*</b>	200 ul	500 ul	
<b>Tissue</b>	≤ 30 mg	500 ul	Homogenization of the tissue in Binding Buffer.
<b>Tissue</b>	≤ 2x10 <sup>6</sup> cells	500ul	Collect up to 2 x 10 <sup>2</sup> cells. Resuspend the pellet in Binding Buffer.
<b>Milk</b>	100 ul	600**ul	** Dilution of 100 ul of milk in 600 ul of working solution (500 ul PV1 supplemented with 100 ul of sample preparation buffer SPM)

\*EDTA-blood, serum, amniotic fluid, CSF, urine, water, etc.

## 6.2 Procedure

The following procedures are for the preparation of 200 µl nucleic acids of sample volume. If larger sample volumes (up to 300 µl) or other sample matrices are used, refer to Table 3 for appropriate buffer volumes.

Samples containing precipitates must be centrifuged before purification!

Store the eluted RNA at ≤ 65 ° C and the eluted DNA at ≤ -18 ° C for further analysis.

## 6.3. Magnetic beads management

### 6.3.1. Magnetic beads distribution

A homogeneous distribution of the magnetic beads in the individual wells of the separation plate is essential for high consistency. Therefore, before distributing the marbles, make sure the marbles are completely homogenized. Shake the bottle well or put it in a vortex. The premixing magnetic beads with the binding buffer allow an easier homogeneous distribution of the marbles in the individual wells of the separation plate. During automation, it is advisable to perform a premix phase before sucking the marbles / collection buffer mixture from the tank to keep them resuspended.

### 6.3.2. Magnetic separation time

The attraction of the magnetic beads to the magnetic pins depends on the magnetic force of the magnetic pins, the selected separation plate, the distance of the separation plate from the magnetic pins and the volume to be processed. The individual times for the complete attraction of the balls to the magnetic pins must be checked and adjusted on each system. It

## EXMB003 INTEGRA RNA/DNA EXTRACTION KIT

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

is recommended to use the separation plates or tubes specified by the magnetic separator supplier.

### 6.3.3. Beads washing

The washing of the marbles can be obtained by shaking or mixing. In contrast to the pipetted mixing up and down, mixing by stirrer or magnetic mixing allows simultaneous mixing of all samples. This reduces the time and the number of tips needed for preparation. Resuspension by pipetting is however more effective than mixing from an agitator or a magnetic mix.

### 6.4. Genomic DNA/RNA extraction

This protocol must be followed if necessary to simultaneously extract genomic DNA, viral DNA and RNA, bacterial DNA or DNA and RNA. Simultaneous isolation of DNA (bacterial or viral) and viral RNA is recommended if eluates are to be used, e.g. for real-time multiplex PCR (RT-) for the detection of both pathogens with DNA genome and RNA virus. For the (simultaneous) isolation of DNA and / or viral RNA from liquid sample material, follow the protocol below for the isolation of viral DNA and / or RNA from other types of samples or complex sample matrices, making refer to the overview in Table 3 or contact us.

Before starting, prepare a working solution of the binding buffer (PV1) integrated with reconstituted poly A (PA) for at least one sample (N) more than necessary to compensate the foaming of the buffer.

**Table 4 – Preparation of working solutions of DNA and RNA.**

Sample needed volume	Working solution of the mastermix
500 ul Binding Buffer (PV1)	500 ul x (N+1)
4 ul Poly A (PA)	4 ul x (N+1)

**Preheat (60 ° C) the Elution Buffer (P4) can increase the result in DNA.**

#### 6.4.1. Manual protocol

This protocol is for use and serves as a guideline for adapting the kit to robotic instruments.

##### Step 1

- Add a working solution of 500 µl, freshly prepared in a 2.0 ml microcentrifuge tube without nuclease.
- Add the 200 µl sample to the microcentrifuge tube.
- Mix immediately.
- Incubate for 60 minutes at 56-60 ° C.
- After lysis incubation, centrifuge for 5 seconds at max. g to collect any sample from the lysis tube and transfer each lysate into the wells of a Square-well Block.

##### Step 2

- Add 20 µl of INTEGRA (MB) magnetic beads to the lysate.
- Stir immediately.
- Incubate for 10 minutes at room temperature with stirring (optional mix by pipetting up and down).
- Separate the magnetic beads from the side of the wells, positioning the Square-well Block on a magnetic separator. Wait at least 30 seconds until all the marbles have been attracted

## EXMB003 INTEGRA RNA/DNA EXTRACTION KIT

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

to the magnets. Remove and discard the supernatant by pipetting. Do not touch the attracted marbles by aspirating the supernatant.

### **Step 3**

- Remove the Square-well Block from the magnetic separator.
- Add 500 µl of inhibitor removal buffer (P2) and resuspend the beads by shaking (optional mixing by pipetting up and down) until the beads are completely resuspended (at least 30 seconds).
- Separate the magnetic beads from the side of the wells by inserting the Square-well Block of the magnetic separator. Wait at least 30 seconds until all the marbles have been attracted to the magnets. Remove and discard the supernatant by pipetting. Do not disturb the attracted marbles by aspirating the supernatant.

### **Step 4**

- Remove the Square-well Block from the magnetic separator.
  - Add a 450 µl Wash buffer (P3) and resuspend the marbles by shaking (optional mixing by pipetting upwards and the city) until the beads are completely repositioned (at least 30 seconds)
  - Separate the magnetic beads from the side of the wells by inserting the Square-well Block of the magnetic separator. Wait at least 30 seconds until all the marbles have been attracted to the magnets. Remove and discard the supernatant by pipetting. Do not disturb the attracted marbles by aspirating the supernatant.
- \* Optional: step 4 can be repeated to increase the degree of purity.

### **Step 5**

- Dry the magnetic ball solution for 5-10 minutes at room temperature.

### **Step 6**

- Remove the Square-well Block from the magnetic separator.
  - Add 70-100 µl (optional: preheated [56-60 ° C]) Elution buffer (P4). Incubate for 10 minutes at room temperature with stirring.
  - Separate the magnetic beads from the side of the wells by inserting the Square-well Block of the magnetic separator. Wait at least 30 seconds until all the marbles have been attracted to the magnets.
  - The supernatant contains purified nucleic acid.
  - Transfer the supernatant to elution plates.
- \* For complex sample types, such as whole blood or stool samples, a further wash step with the Wash buffer (P3) can increase the purity of the nucleic acid.

## **6.4.2. Automatic protocol on KingFisher™ Flex Magnetic Particle Processor**

### **Step 1 – Sample preparation Part I and Lysis**

- **Add 500 µl** of freshly prepared work solution to each well of a 96-well plate.
- **Add 200 µL** of sample.
- **Mix** immediately.
- **Incubate** for 60 minutes at 56-60°C.

### **Step 2 – Washing plates preparation**

- **Add 500 µl** of Inhibitor removal buffer (P2) to each well of a 96 plate deep-well.

## EXMB003 INTEGRA RNA/DNA EXTRACTION KIT

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

- Add 450 µl of Wash buffer (P3) to each well of the deep-well plate.
- Add 450 µl of Wash buffer (P3) to each well of a second deep-well plate.

### **Step 3 – Elution plate preparation**

- Add 70-100 µl of Elution buffer (P4) of each well of a plate with 96 deep-well.

### **Step 4 – sample preparation Part II/binding**

- Add 20 µl of magnetic particles *INTEGRA* to the lysate.

### **Step 6 – Purification protocol execution on the equipment**

- Insert the plates as indicated by the KingFisher™ tool.
- the method starts with a combined lysis and binding step after setting the last plate on the instrument.

### **Step 7 – removal of eluted RNA/DNA**

- The instrument stops after the last elution step. Following the instructions of the instrument, download the used plates.
- Purified RNA / DNA can now be used for detection.

For the KingFisher™ Flex Magnetic Particle Processors settings, use the guide that can be supplied separately from IMMUNOSPARK.

## **6.5. DNA extraction from raw milk**

Needed additional materials:

- Sample Preparation Buffer SPM

Before starting, prepare a working solution of the release buffer (PV1) integrated with reconstituted polymer (PA) and SPM sample preparation buffer for at least one sample (N) more than necessary to compensate the foaming of the buffer.

**Table 5 – Working solution preparation for extraction of DNA and RNA from milk.**

Needed volume per sample	Mastermix working solution
500 ul Binding Buffer (PV1)	500 ul x (N+1)
4 ul Poly A (PA)	4 ul x (N+1)
100 ul Samples Preparation Buffer	100 ul (N+1)

### **Step 1**

- Add a working solution of 600 µl, freshly prepared in a 2.0 ml microcentrifuge tube without nuclease.
- Add 100 µl of sample to each sample.
- Mix immediately.
- Digestion for 30 minutes at 60 ° C.
- After lysis incubation, centrifuge for 5 seconds at max. speed to collect any sample from the lysis tube lids.

**Step 2** and **following** see the instruction for the extraction of genomic DNA/RNA.

## **7. TECHNICAL INFORMATION**

### **7.1 Quality control**

According to the ISO procedures of Immunospark, each lot of *INTEGRA* is tested on the basis of pre-established specifications to guarantee a constant product quality.

## EXMB003 INTEGRA RNA/DNA EXTRACTION KIT

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

The dilution series of MS2 RNA, Lambda-Phage DNA and human genomic DNA are applied to INTEGRA marbles, washed and eluted according to the kit protocol. 4 µl of the eluate is analyzed in real time (RT-) PCR. Recovery of at least  $2 \times 10^5$  RNA or DNA molecules for 200 µl of sample is guaranteed.

For the validation of the viral RNA of INTEGRA, genomic DNA, viral DNA and bacterial DNA have been isolated from a wide range of sample matrices such as blood, tissues, feces, ticks, milk, buccal swab, etc.

The eluates were used as real-time RT-PCR templates and produced highly specific PCR products with good yields.

### 7.2 Preparation time

The required preparation time always depends on the number of samples to be prepared.

**Table 6 – Preparation time.**

	DNA/RNA	Viral RNA
<b>Total time</b>	Aprox. 100 minutes	Aprox. 40 minutes
<b>Manual time</b>	Less than 10 minutes	Less than 10 minutes

## 8. TROUBLESHOOTING

The following troubleshooting guide is included to help you with any problems that may arise during the isolation of the nucleic acid from different types of sampling material. Especially when working with complex sample matrices such as adipose tissue, whole blood or highly contaminated samples, sample preparation can be crucial. For protocols on sample materials not covered in this manual or for further questions about nucleic acid isolation, please contact us.



<b>Low quality of nucleic acids</b>	
Insufficient Elution Buffer volume.	The pellet of beads must be completely covered with Elution Buffer.
Insufficient Performance of the buffer during elution	Remove buffer residues during separation steps completely. Residues reduce efficiency of the following washing and elution steps.
Dried particles	Do not allow the marbles to dry out as this may result in decreased efficiency of elution.
Accidental suction of separate magnetic beads	Do not disturb the marbles attracted during aspiration of the supernatant, especially when the magnetic pellet is not visible in the lysate.
Aspiration and loss of magnetic beads	Time for magnetic separation too short or suction speed too high.
Insufficient washing	Use only the appropriate separator combinations and of the plates. Ensure that the beads are resuspended during the washing. If the agitation is not sufficient to resuspend completely, mix repeatedly.
Carry over of ethanol from Wash buffer	Be sure to remove all the ethanol Wash buffer from the final wash, because it interferes.
Evaporation of ethanol from Wash buffer	Close the bottles of the buffer well, avoiding the evaporation of the ethanol. Do not reuse the pads.
Magnetic separation time too short	Increase the separation time to allow the marbles to be completely attracted by the magnetic pins before drawing any liquid from well.

**EXMB003 INTEGRA RNA/DNA EXTRACTION KIT**

Info: [sales@immunospark.com](mailto:sales@immunospark.com)

Excessive suction speed (Elution)	The high suction speed during the elution step can cause the suction of the balls. Reduce the suction speed for the elution step.
Kit improperly stored (T)	Store the kit at +18 a + 25°C at the arrival.
Reagents exposed to external agents	Store all buffers at +18 to + 25 ° C. Close all reagent bottles well after each use to preserve pH and stability and to prevent contamination. Aliquot Poly A (PA) after reconstitution and store aliquots at ≤ -18 ° C.
2-Propanol not added to PV1	Add 2-propanol to the buffer before using it. Mix the buffer well and store at +18 at 25 ° C. Always mark the buffer bottle to indicate whether or not 2-propanol was added.
Ethanol not added to P2 e/o P3	Add absolute ethanol to the buffers before use. After adding ethanol, mix the buffers well and store at +18 to 25 ° C. Always mark the buffer bottles to indicate whether or not ethanol has been added.
Reagents and samples not completely mixed	Always mix the sample tube well after adding each reagent.
Impurities not completely removed	Perform a second wash step with the Wash buffer (P3) to completely remove salts, proteins and other residual impurities from the bound nucleic acid.

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