



sFlt-1
PlGF

Pre-eclampsia management with biomarkers

Improved pre-eclampsia diagnosis and prognosis of adverse outcome with PlGF and sFlt-1 after week 20 of gestation


Biomarkers for pre-eclampsia management

Improving the diagnostic tools for pre-eclampsia evaluation

Pre-eclampsia is a progressive, pregnancy-related disorder with severe complications for mother and child. A timely diagnosis is needed in order to prevent maternal and fetal morbidity or mortality. In the absence of a specific therapy other than delivery the main objective of frequent patient monitoring is to detect deterioration of a patient's condition and to counteract maternal and fetal risk.



Diagnostic standard for pre-eclampsia ...



Blood pressure



Proteinuria



10% of pregnant women show unspecific signs and symptoms of pre-eclampsia



Only **one fifth of them** is actually developing pre-eclampsia¹


The “gold standard” for pre-eclampsia diagnosis – assessment of blood pressure and proteinuria – offers only a **poor sensitivity and specificity** with regards to origin of disease and prediction of maternal and perinatal outcome.^{2,3}

Serum sFlt-1 and PIGF determination adds significant clinical benefit to standard procedures

Determination of the biomarkers sFlt-1 (soluble FMS-like Tyrosine Kinase) and PIGF (Placental Growth Factor) in maternal blood have shown to significantly improve risk stratification among women presenting for pre-eclampsia evaluation.

The high sensitivity assays Thermo Scientific™ B·R·A·H·M·S™ sFlt-1 KRYPTOR™ and Thermo Scientific B·R·A·H·M·S PIGF plus KRYPTOR detect serum levels of both biomarkers reliably throughout pregnancy and thus improve pre-eclampsia management.

... significantly improved by sFlt-1 and PIGF serum measurement



Clearly **differentiating** between pre-eclampsia and other forms of hypertensive disorders^{4,5}



Reliable **prognosis** of adverse outcome in women with suspected pre-eclampsia^{5,6}



Offering **potential savings** in hospital costs and resource use⁷

Measuring sFlt-1 and PIGF starting in mid pregnancy in women with suspected pre-eclampsia significantly improves the current evaluation of patients – for a better patient management and improved care.

Pre-eclampsia diagnosis and prognosis of adverse outcome

The added value of sFlt-1 and PIGF


Improved diagnosis of pre-eclampsia with sFlt-1/PIGF ratio

Studies proved the additional benefit of the sFlt-1/PIGF ratio in diagnosing pre-eclampsia:

- In women presenting with hypertension, the sFlt-1/PIGF ratio is able to **distinguish between those who will develop pre-eclampsia and those with chronic or gestational hypertension**. Women with pre-eclampsia have a significantly higher sFlt-1/PIGF ratio than women with other hypertensive disorders or controls.^{4,5}
- The addition of sFlt-1/PIGF ratio to Doppler ultrasound measurement **improves the sensitivity and specificity** in diagnosing pre-eclampsia compared to the Doppler measurement alone.⁶

- Measurement of sFlt-1 and PIGF levels in maternal serum, starting in mid pregnancy, can **confirm pre-eclampsia diagnosis**, with the sFlt-1/PIGF ratio having a superior diagnostic ability compared to either of the biomarkers alone.^{8,9}

Therefore, the sFlt-1/PIGF ratio is a valuable additional tool for confirming or excluding the diagnosis of pre-eclampsia.



PIGF and sFlt-1 were measured on KRYPTOR in parallel on samples from pregnant women with normal pregnancy outcome and patients with pre-eclampsia. At a cut-off of 85 for the sFlt-1/PIGF ratio, the sensitivity was calculated at 95% and the specificity at 84% for diagnosing pre-eclampsia. Latest studies show identical clinical performance and high accuracy in diagnosing pre-eclampsia when applying recently published cut-offs by using the sFlt-1/PIGF ratio on KRYPTOR.^{6,10,11}

The higher the sensitivity of a test the more women with pre-eclampsia are identified correctly and can be advised for closer monitoring.

Figure 1 Improved pre-eclampsia diagnosis with sFlt-1/PIGF ratio¹²



Prognosis of adverse outcome with sFit-1/PIGF ratio

Recent studies showed that **women with any subsequent adverse outcome** in addition to hypertension had a significantly higher sFit-1/PIGF ratio than those women without, especially when presenting before week 34 (Figure 2).^{5,6}




Figure 2 Prediction of adverse outcome with sFit-1/PIGF ratio in women presenting < 34 weeks' gestation⁶

Women who needed to be delivered within the next 2 weeks after presentation had a significantly higher sFit-1/PIGF ratio than women who could continue with their pregnancy (Figure 3).^{5,6}




Figure 3 Prediction of duration of pregnancy with sFit-1/PIGF ratio in women presenting < 34 weeks' gestation⁶


The sFit-1/PIGF ratio is also a potent predictor for subsequent maternal and fetal adverse outcome in women already diagnosed with pre-eclampsia and can support clinical decisions.

The role of angiogenic factors

Biomarker levels correlate with severity of disease

sFlt-1 and PlGF are counterparts

Although the cause of pre-eclampsia remains unclear, it is likely that the syndrome is initiated by an imbalance of angiogenic factors secreted by the placenta that induce endothelial dysfunction.



sFlt-1 is a truncated form of the VEGF receptor Flt-1, circulating freely in the blood. sFlt-1 is produced in the placenta and secreted into the bloodstream, where it binds VEGF and PlGF with high affinity and therefore neutralizes their effects.⁸



PlGF belongs to the Vascular Endothelial Growth Factors (VEGF) family, promoting proliferation and survival of endothelial cells and inducing vascular permeability.¹³

••▶ Signal transduction (healthy)

▶ Signal transduction inhibited

Figure 4 sFlt-1 acts as potent antagonist of PlGF and VEGF by adhering to the receptor-binding domains, thus preventing interaction with endothelial receptors and inducing endothelial dysfunction





Angiogenic factors during pregnancy

Normal pregnancy

During pregnancy, sFlt-1 levels are stable until weeks 20–24, when they rise steadily until delivery. In contrast, PIGF levels increase progressively in first and second trimester and decrease towards term.¹³

Pre-eclamptic pregnancy

In women with pre-eclampsia, sFlt-1 levels are significantly increased while concentrations of circulating free PIGF are significantly decreased.^{13,14} In contrast to PIGF where the difference between healthy and pre-eclamptic pregnancies is measurable throughout pregnancy, sFlt-1 levels only start to separate after week 20.

- Controls
- Women who had pre-eclampsia >5 weeks later
- Women who later had pre-eclampsia
- Women with clinical pre-eclampsia

Figure 5 Mean sFlt-1 and PIGF concentrations of healthy women and those women who later developed pre-eclampsia¹³

Measuring maternal serum concentrations of sFlt-1 and PIGF can differentiate healthy women from women with pre-eclampsia.^{9,15} Changes in sFlt-1 and PIGF levels also reflect the severity of the disease: early-onset pre-eclampsia is associated with greater changes compared to late-onset pre-eclampsia.¹⁶

Pre-eclampsia management throughout pregnancy

Improving the outcomes for mother and child

PIGF and PAPP-A: First trimester screening for timely intervention



Combined screening for pre-eclampsia in weeks 11–13+6 can reliably identify women at risk for developing pre-eclampsia.

Combined first trimester screening includes

- serum PIGF and PAPP-A measurement,
- determination of mean arterial pressure (MAP), and
- Uterine Artery Pulsatility Index (UtA-PI)

resulting in a detection rate of >90% for a fixed false positive rate of 5%.¹⁷

An early identification of high-risk women allows for preventive measures and intensified monitoring. Administering low-dose aspirin (<150 mg/day) to high-risk women before 16 weeks of gestation can significantly reduce the incidence of pre-eclampsia by 50%–90%.^{18,19}



Facts on pre-eclampsia

- Multisystem, life-threatening pregnancy-related disorder
- A main reason for maternal and fetal morbidity and mortality^{20,21}
- **Incidence:** 2-8% of pregnancies²⁰
- **Definition:** New onset hypertension and proteinuria >20 weeks of gestation in previously normotensive women²²
- **HELLP** (Hemolysis, Elevated Liver enzymes, Low Platelets): Severe pre-eclampsia variant occurring in ~20% of symptomatic women; defined by additional affection of liver and coagulation system²³
- **Eclampsia:** Final stage of disease, associated with severe tonic-clonic seizures and coma as well as brain injury, cerebral edema and stroke²³

sFlt-1/PIGF ratio: Improved diagnosis and prognosis of adverse outcome

First symptoms of pre-eclampsia (hypertension, proteinuria) are observed after 20 weeks of gestation.²³

Diagnosis of pre-eclampsia is difficult, as pre-eclampsia can be confused with other diseases such as pregnancy-induced hypertension.

By adding sFlt-1/PIGF ratio to the current diagnostic standard, the **diagnosis of pre-eclampsia** in a symptomatic woman can be confirmed or excluded.^{2,10,11}

In women with diagnosed pre-eclampsia, the sFlt-1/PIGF ratio is a potent **predictor of subsequent maternal and fetal adverse outcome** and can be useful for further patient management.^{5,6}




Figure 6 First clinical symptoms of pre-eclampsia are observed >20 weeks of gestation. The gestational age at onset correlates with the severity of maternal and fetal consequences.²



Complete pre-eclampsia portfolio

From safe screening to improved diagnosis with B·R·A·H·M·S sFlt-1 and PIGF plus

Thermo Scientific B·R·A·H·M·S sFlt-1 KRYPTOR

Automated immunofluorescent assay for the quantitative determination of the concentration of sFlt-1 (soluble FMS-like Tyrosine Kinase 1, also known as VEGF receptor-1) in human serum.

- 75 determinations per kit
- 9 min incubation time
- Monoparametric control kit, 3 levels
- Wide measuring range: 22–90 000 pg/mL
- Excellent precision

With the lower and upper detection limits of 22 and 90 000 pg/mL B·R·A·H·M·S sFlt-1 KRYPTOR provides the measuring range needed for a **reliable detection of clinical sFlt-1 values throughout pregnancy.**

Thermo Scientific B·R·A·H·M·S PIGF plus KRYPTOR

Automated immunofluorescent assay for the quantitative determination of the concentration of PIGF (Placental Growth Factor) in human serum. The assay is specific for the measurement of human free PIGF-1.

- 75 determinations per kit
- 29 min incubation time
- Monoparametric control kit, 3 levels
- Wide measuring range: 3.6–7 000 pg/mL
- Excellent precision

With a detection limit of 3.6 pg/mL and an upper limit of 7 000 pg/mL B·R·A·H·M·S PIGF plus KRYPTOR provides the high sensitivity needed for **measuring PIGF levels in first trimester** as well as a wide measuring range to **reliably measure clinical values throughout pregnancy.**




sFlt-1
PIGF

Exceptionally precise, fast and easy

B·R·A·H·M·S KRYPTOR analyzers

- All KRYPTOR platforms FMF approved
- In routine use by FMF since 1999
- Excellent precision and proven median stability



Thermo Scientific
B·R·A·H·M·S KRYPTOR GOLD

Thermo Scientific
B·R·A·H·M·S KRYPTOR compact PLUS



Your **BENEFITS** by adding the sFit-1/PIGF ratio into clinical routine

- ▶ **Improved clinical accuracy** in diagnosing pre-eclampsia in symptomatic patients
- ▶ **Potent prognostic tool** for subsequent adverse pregnancy outcome

Your **ACCESS** to our interactive e-detail

Get more information on pre-eclampsia management throughout pregnancy:



<http://prenatal.world-of-biomarkers.com>

Pin code: **ratio01**



Thermo Scientific B·R·A·H·M·S Biomarkers Prenatal Screening Portfolio on B·R·A·H·M·S KRYPTOR Systems

B·R·A·H·M·S AFP KRYPTOR	Art. no.: 816.075
B·R·A·H·M·S Free βhCG KRYPTOR	Art. no.: 809.075
B·R·A·H·M·S hCG+β KRYPTOR	Art. no.: 841.050
B·R·A·H·M·S Inhibin A KRYPTOR*	Art. no.: 850.075
B·R·A·H·M·S PAPP-A KRYPTOR	Art. no.: 866.075
B·R·A·H·M·S PIGF plus KRYPTOR**	Art. no.: 859.075
B·R·A·H·M·S sFit-1 KRYPTOR**	Art. no.: 845.075
B·R·A·H·M·S uE3 KRYPTOR***	Art. no.: 803.075
B·R·A·H·M·S Fast Screen pre I plus Software	Art. no.: 105750

* Available on B·R·A·H·M·S KRYPTOR GOLD

** Available on B·R·A·H·M·S KRYPTOR compact PLUS and B·R·A·H·M·S KRYPTOR GOLD

*** Available on B·R·A·H·M·S KRYPTOR, B·R·A·H·M·S KRYPTOR compact PLUS and B·R·A·H·M·S KRYPTOR GOLD

References

- Milne et al. *BMJ* 2009; 339: b3129
- Verlohren et al. *Clin Sci* 2012; 122 (2): 43-52
- Zhang et al. *Obstet Gynecol* 2001; 97: 261-7
- Verlohren et al. *Am J Obstet Gynecol* 2012; 206, 58.e1-8
- Rana et al. *Circulation* 2012; 125 (7): 911-9
- Salahuddin et al. *Hypert Preg* 2016; 35 (3): 330-45
- Schnettler et al. *BJOG* 2013; 120: 1224-32
- Hagmann et al. *Clin Chem* 2012; 58 (5), 837-45
- Verlohren et al. *Am J Obstet Gynecol* 2010; 202: 161.e1-11
- Andersen et al. *J Am Soc Hypertens* 2015; 9 (2): 86-96
- van Helden et al. *Clin Biochem* 2015; 48 (16-17): 1113-9
- B·R·A·H·M·S sFit-1 KRYPTOR Instructions for use, 2013
- Levine et al. *N Engl J Med* 2004; 350: 672-83
- De Vivo et al. *Acta Obstet et Gynecol* 2008; 87: 837-42
- Thadhani et al. *J Clin Endo* 2004; 89 (2): 770-5
- Romero et al. *J Matern fetal Neonatal Med* 2008; 21: 9-23
- Akolekar et al. *Prenat Diagn* 2011; 31: 66-74
- Bujold et al. *Obstet Gynecol* 2010; 116: 402-14
- Park et al. *Ultrasound Obstet Gynecol* 2015; 46: 419-423
- WHO. *World Health Report* 2005: 63
- Duley et al. *Semin Perinatol* 2009; 33 (3): 130-7
- Definition of the American College of Obstetrics and Gynecology
- Powe et al. *Circulation* 2011; 123: 2856-69
- Poon et al. *Hypertension* 2009; 53: 812-8

Clinical Diagnostics	Thermo Fisher Scientific	+49 (0)3302 883 0
	B·R·A·H·M·S GmbH	+49 (0)3302 883 100 fax
	Neuendorfstr. 25	info.brahms@thermofisher.com
	16761 Hennigsdorf Germany	www.thermoscientific.com/brahms

Find out more at thermoscientific.com/brahms

Products are CE marked but not 510(k)-cleared and not available for sale in the U.S. Availability of products in each country depends on local regulatory marketing authorization status.

© 2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. KRYPTOR is a trademark of Cisbio Bioassays, licensed for use by B·R·A·H·M·S, a part of Thermo Fisher Scientific.

107304.4

ThermoFisher
SCIENTIFIC