

KRYPTOR : FERTILITY



	PROLACTIN	hCG+β
Available on KRYPTOR compact plus	NO	OUI
Sample volume (µL)	26	26
Nb of tests per kit	50	50
Reagent kit	Ready-to-use	Freeze-dried (Cr)
Incubation time (min)	19	14
Automatic dilution factor	5, 10, 50 et 300	10 100 1000
Direct measurement range	0 - 4600 µL/mL	2 - 2000 UI/L
Normal values	Men 5e - 95e percentile 58 - 254 Non menopausal women : 80 - 435 Menopausal women : 68 - 318 Children Non pubescent: 43 - 318	Women * Non pregnant : 100% < 5 UI/L * Menopausal : 98 % < 5 UI/L men 100 % < 5 UI/L
Detection limit	5 µUL/mL	2 UI/L
On-board stability (days)	29	29
Calibrators	1 Freeze-dried	1
Calibration frequency (days)	15	15
International standards	3 ^{ème} IS OMS 84/500	3 ^{ème} IS OMS 75/537