

KRYPTOR : INFECTION



	PROCALCITONINE (PCT)
Available on KRYPTOR compact plus	YES
Sample volume (µL)	50
Nb of tests per kit	50
Reagent kit	Freeze-dried
Incubation time (min)	19
Automatic dilution factor	20 100
Direct measurement range	0.02 – 50 µg/L
Normal values	<0.06 µg/L
Detection limit	0.02 µg/L
On-board stability (days)	29
Calibrators	1 Freeze-dried
Calibration frequency (days)	15
International standards	-