


KRYPTOR : OTHERS PRODUCTS



	FERRITIN	OSTEOCALCIN	CRP US
Available on KRYPTOR compact plus	YES	YES	NO
Sample volume (µL)	26	8	10 µl
Nb of tests per kit	50	50	75
Reagent kit	Ready-to-use	Ready-to-use	Ready-to-use
Incubation time (min)	19	14	9
Automatic dilution factor	10 100 500	10	10
Direct measurement range	0.8 - 1200 ng/mL	0.4 - 270 ng/mL	0.06 – 200 µg/mL
Normal values	95 th percentile Men 274 Non menauposal women 121	Extreme values Men 5.1-25.2 women 6.1-43.5	-
Detection limit	0.8 ng/mL	0.4 ng/mL	<0.06 µg/mL
On-board stability (days)	29	29	29
Calibrators	1	1 Freeze-dried	2 (2 calibrators for DF measurement and 1 calibrator for ratio measurement)
Calibration frequency (days)	15	15	15
International standards	1 ng CIS = 1 UI 2 nd IRP OMS 80/578 1 ng CIS = 1.05 UI 1 st IRP OMS 80/602		0.06 µg/mL